

TECHNICAL INFORMATION

Cabinet Mounted Controls



Panel Mount Auto Control
Model #29-355

STANDARD FEATURES

- 3 Handpiece control
- Automatic control
- Individual handpiece coolant controls
- Handpiece purge system
- Asepsis handpiece tubings
- Master on/off control
- Wet/dry foot control
- 3-way air/water syringe
- Individual handpiece pressure adjustment



Micro II Manual
Model #28-101

STANDARD FEATURES

- 2 Handpiece control
- Asepsis handpiece tubings
- Standard foot control
- 3-way air/water syringe
- Individual handpiece pressure adjustment



Micro II Auto
Model #28-115

STANDARD FEATURES

- 2 Handpiece control
- Automatic Control
- Asepsis handpiece tubings
- Standard foot control
- 3-way air/water syringe
- Individual handpiece pressure adjustment