

TECHNICAL INFORMATION

Cabinet Mounted Systems



*12 O'clock
Dual Purpose System
AS-1 Unit w/Assistant's Arm
Model #70-900*

STANDARD FEATURES

- 3 Handpiece control
- Automatic control
- Individual handpiece coolant controls
- Handpiece purge system
- Asepsis handpiece tubings
- Master on/off control
- Reversible from right to left
- Junction box w/ master valves, filters & regulators
- Wet/dry foot control
- 3-way air/water syringe
- Individual handpiece pressure adjustment
- Horizontal travel arm with dual pivot
- Solids collector
- HVE
- Saliva ejector
- Assistant's air/water syringe
- Assistant's 3 Position Holder Bar



*12 O'clock
Doctor's System
AS-1 Unit only
Model #70-902*

STANDARD FEATURES

- 3 Handpiece control
- Automatic control
- Individual handpiece coolant controls
- Handpiece purge system
- Asepsis handpiece tubings
- Master on/off control
- Reversible from right to left
- Junction box w/ master valves, filters & regulators
- Wet/dry foot control
- 3-way air/water syringe
- Individual handpiece pressure adjustment
- Horizontal travel arm with dual pivot



*12 O'clock
Assistant's System
Model #70-904*

STANDARD FEATURES

- Reversible from right to left
- 3-way air/water syringe
- Horizontal travel arm with dual pivot
- Solids collector
- HVE
- Saliva ejector
- Assistant's 3 Position Holder Bar

OPTIONAL ACCESSORIES

- Standard handpiece illumination system
- 6-Pin handpiece illumination system
- 4th handpiece station
- Junction box w/ master valves, filters & regulators
- Clean water system
- Additional HVE