

TECHNICAL INFORMATION

Cuspidor & Assistant's Arm



Porcelain Cuspidor
Model #69-240

STANDARD FEATURES

- Porcelain cuspidor
- Adjustable timed bowl rinse
- Push button activation for cup fill and bowl rinse
- Asepsis design for easy cleaning
- Water outlet
- Selector switch for clean/city water
- Solids collector
- HVE
- Saliva ejector
- Clean water system
- Standard 5ft umbilical



Telescoping Arm
Model #69-201

STANDARD FEATURES:

- Asepsis design for easy cleaning
- Water outlet
- Selector switch for clean/city water
- 3-Way air/water syringe
- Solids collector
- HVE
- Saliva ejector
- Clean water system
- Telescoping arm with 4 position holder bar
- Additional HVE
- Standard 5ft umbilical



AS Style Cuspidor
Model #69-202

STANDARD FEATURES:

- Removable glass bowl
- Adjustable timed bowl rinse
- Push button activation for cup fill and bowl rinse
- Asepsis design for easy cleaning
- Water outlet
- Selector switch for clean/city water
- Solids collector
- HVE
- Saliva ejector
- Clean water system
- Standard 5ft umbilical

OPTIONAL ACCESSORIES

- 3-Way air/water syringe for cuspidor
- Junction box w/master valves, filters & regulators
- Telescoping Assistant's Arm
- Additional umbilical
- Chair bracket
- Mounting post (8-1/4")