

# TECHNICAL INFORMATION

## Engle 320 with ComfortPlus™



"ComfortPlus™" chair shown with optional Ultra Leather upholstery with "TwinStitch™"

New "ComfortPlus™" provides Synchronized Anatomically Correct Positioning to maximize the patient's comfort.

### ENGLE 300 HYDRAULIC CHAIR

#### Standard Features

- Thin back with lumbar support
- Anatomically correct back pivot point
- Adjustable, quick release double articulating headrest with contour cushion
- Four manual chair movements and five programmable positions are standard along with a programmable exit
- Heavy vinyl scuff guard for toe board
- Sturdy all-metal construction
- 60° chair rotation
- All pivot points have permanently lubricated bearings
- Silent hydraulic lift system with auto rise safety shut-off
- Electronic system 10 year limited warranty

### ENGLE 320 DELIVERY SYSTEM

- Delivery system provides effortless flex-arm with air lock positioning
- Unit mounted touch pad
- Clean water system
- 3 Handpiece positions will accommodate extra lines for ancillary products
- Wet/Dry foot control
- Horizontal travel arm with assistant's instrumentation

### LIGHT SYSTEM

- 3 axis of rotation
- Sharply defined light pattern
- Dual intensity settings
- Ideal color temperature (4500°K)

Note: Chair shown with optional TwinStitch™ with Ultra Leather and internal umbilical.