

TECHNICAL INFORMATION

Hygiene Systems



*Hygiene Unit
with Cuspidor
Model #70-885*

STANDARD FEATURES

- 2 HP control
- Syringe
- Purge system
- Independent coolant controls
- Clean water system
- Counter-balanced flex arm with air brake
- Swivel tray

UNITS WITH CUSPIDOR INCLUDE

- Solids collector
- Saliva ejector
- Timed bowl rinse
- Cup filler
- Clean water system

OPTIONAL ACCESSORIES

- Handpiece illumination system

Standard Features

- 2 HP control
- Syringe
- Purge system
- Independent coolant controls
- Clean water system
- Counter-balanced flex arm with air brake
- Swivel tray

TELESCOPING ARM

- 2 HVE
- Saliva ejector
- Air/water syringe
- Water outlet
- Solids collector

OPTIONAL ACCESSORIES

- Handpiece illumination system



*Hygiene Unit with
Assistant's Arm
Model #70-886*