

TECHNICAL INFORMATION

Mobile Cart Systems



*Dual Purpose Cart
Model #50-135*

STANDARD FEATURES - MODEL 50-135

- 3 handpiece automatic control
- Individual handpiece coolant controls
- Handpiece purge system
- Asepsis handpiece tubings
- 3-way air/water syringe
- Wet/Dry foot control
- Master on-off switch
- Clean water system
- Duplex electrical water outlet
- Solids collector
- HVE
- Saliva ejector
- Assistant's air/water syringe
- Water quick disconnect
- Junction box w/master valves, filters & regulators
- Storage bins



*Assistant's Cart
Model #50-145*

STANDARD FEATURES - MODEL 50-145

- Master on-off switch
- Clean water system
- Duplex electrical water outlet
- Solids collector
- HVE
- Saliva ejector
- Assistant's air/water syringe
- Water quick disconnect
- Junction box w/master valves, filters & regulators
- Storage bins



*Flip Top Doctor's Cart
Model #50-140*

STANDARD FEATURES - MODEL 50-140

- 3 handpiece automatic control
- Individual handpiece coolant controls
- Handpiece purge system
- Asepsis handpiece tubings
- 3-way air/water syringe
- Wet/Dry foot control
- Master on-off switch
- Clean water system
- Duplex electrical water outlet
- Water quick disconnect
- Junction box w/master valves, filters & regulators



*Flip Top Dual Purpose Cart
Model #50-150*

STANDARD FEATURES - MODEL 50-150

- 3 handpiece automatic control
- Individual handpiece coolant controls
- Handpiece purge system
- Asepsis handpiece tubings
- 3-way air/water syringe
- Wet/Dry foot control
- Master on-off switch
- Clean water system
- Duplex electrical water outlet
- Solids collector
- HVE
- Saliva ejector
- Water quick disconnect
- Junction box w/master valves, filters & regulators

OPTIONAL ACCESSORIES

- Standard handpiece illumination system
- 6-Pin handpiece illumination system
- Air quick disconnect
- 4th handpiece station (AS-2)
- 4th handpiece station (AS-1)
- Swivel tray (for 50-135 & 50-145)
- Additional HVE