

TECHNICAL INFORMATION

Post/Wall Mounted Systems



AS-1
Model #70-875

STANDARD FEATURES

- AS-1 control system
- Three handpiece control system
- Individual coolant controls
- Master on/off controls
- Flex arm with air lock
- Adjustable swivel tray
- 3-way air/water syringe
- Asepsis handpiece tubings
- Handpiece purge system
- Wet/Dry foot control
- Junction box w/ master valves, filters & regulators



AS-2
Model #70-882

STANDARD FEATURES

- AS-2 control system
- Three handpiece control system
- Individual coolant controls
- Master on/off controls
- Flex arm with air lock
- Adjustable swivel tray
- 3-way air/water syringe
- Asepsis handpiece tubings
- Handpiece purge system
- Wet/Dry foot control
- Junction box w/ master valves, filters & regulators



AS-1
Model #70-876

STANDARD FEATURES

- AS-1 control system
- Three handpiece control system
- Individual coolant controls
- Master on/off controls
- Flex arm with air lock
- Adjustable swivel tray
- 3-way air/water syringe
- Asepsis handpiece tubings
- Handpiece purge system
- Wet/Dry foot control
- Junction box w/ master valves, filters & regulators



AS-2
Model #70-883

STANDARD FEATURES

- AS-2 control system
- Three handpiece control system
- Individual coolant controls
- Master on/off controls
- Flex arm with air lock
- Adjustable swivel tray
- 3-way air/water syringe
- Asepsis handpiece tubings
- Handpiece purge system
- Wet/Dry foot control
- Junction box w/ master valves, filters & regulators



12 O'clock
Model #70-892

STANDARD FEATURES

- AS-1 control system
- Three handpiece control system
- Individual coolant controls
- Master on/off controls
- Flex arm with air lock
- Asepsis handpiece tubings
- Handpiece purge system
- Wet/Dry foot control
- Junction box w/ master valves, filters & regulators
- Clean water system
- Saliva collector
- HVE
- Saliva ejector
- Duplex electric outlet
- Water outlet
- Assistant's syringe



12 O'clock
Model #70-893

STANDARD FEATURES

- Master on/off controls
- Flex arm with air lock
- 3-way air/water syringe
- Junction box w/ master valves, filters & regulators
- Clean water system
- Saliva collector
- HVE
- Saliva ejector
- Duplex electric outlet
- Water outlet

OPTIONAL ACCESSORIES

- Standard handpiece illumination system
- 6-Pin handpiece illumination system
- 4th handpiece station (AS-2)
- 4th handpiece station (AS-1)
- Additional assistant's syringe
- Additional HVE
- Storage bins (for 70-892 & 70-893)
- Swivel tray (for 70-892 & 70-893)
- Wall mount board