

TECHNICAL INFORMATION

Standard Doctor Stools



STANDARD DOCTOR STOOL FEATURES

- Five legged base
- Twin wheel casters
- Adjustable back support

STANDARD ASSISTANT STOOL FEATURES

- Five legged base
- Twin wheel casters
- Rotary abdominal support

Deluxe Doctor Stools

DELUXE DOCTOR STOOL FEATURES

- Five legged base
- Twin wheel casters
- Adjustable back support
- Adjustable back tilt
- Adjustable seat tilt

DELUXE ASSISTANT STOOL FEATURES

- Five legged base
- Twin wheel casters
- Rotary abdominal support
- Adjustable footring
- Quick release for entry/exit

