

# TECHNICAL INFORMATION

## Chair Mounted Systems



AS-2 with Utility Center  
Model #70-912

- STANDARD FEATURES**
- AS-2 control system
  - 3 Handpiece automatic system
  - Individual handpiece coolant controls
  - Master on/off switch
  - Flex arm with air lock
  - Adjustable swivel tray
  - 3-way air/water syringe
  - Asepsis handpiece tubings
  - Wet/Dry foot control
  - Utility center
  - Water outlet
  - Clean water system
  - Oversize junction box w/master valves, filters & regulators
  - Chair mounting bracket
  - Light post assembly
  - Selector switch for clean/city water



AS-2 with AS Style Cuspidor  
Model #70-878

- STANDARD FEATURES**
- AS-2 control system
  - 3 Handpiece automatic system
  - Individual handpiece coolant controls
  - Master on/off switch
  - Flex arm with air lock
  - Adjustable swivel tray
  - 3-way air/water syringe
  - Asepsis handpiece tubings
  - Wet/Dry foot control
  - AS Style Cuspidor
  - Utility center
  - HVE
  - Saliva ejector
  - Water outlet
  - Clean water system
  - Oversize junction box w/master valves, filters & regulators
  - Chair mounting bracket
  - Light post assembly
  - Selector switch for clean/city water



Euro Control with Utility Center  
Model #70-914

- STANDARD FEATURES**
- Euro control system
  - 3 Handpiece automatic system
  - Individual handpiece coolant controls
  - Master on/off switch
  - Flex arm with air lock
  - Adjustable swivel tray
  - 3-way air/water syringe
  - Asepsis handpiece tubings
  - Wet/Dry foot control
  - Utility center
  - Water outlet
  - Clean water system
  - Oversize junction box w/master valves, filters & regulators
  - Chair mounting bracket
  - Light post assembly
  - Selector switch for clean/city water

**OPTIONAL ACCESSORIES**

- Standard handpiece illumination system
- 6-Pin handpiece illumination system
- Telescoping Assistant's Arm
- AS-1 Control Head substitution
- 4th handpiece station (excluding Euro system)
- Additional assistant's syringe
- Air quick disconnect
- Additional HVE
- Low voltage switch
- Euro Control Head