

TECHNICAL INFORMATION

Engle 1200 Series Package



ENGLE 1200 SERIES PACKAGE

STANDARD FEATURES

- AS-2 three handpiece control system
- Individual handpiece coolant controls
- Asepsis handpiece tubings
- Wet/Dry foot control
- Handpiece purge system
- Master on/off switch
- Flex arm with air lock
- Adjustable swivel tray
- 3-way air/water syringe
- Porcelain Cuspidor with utility center
- Saliva collector
- HVE
- Saliva ejector
- Clean water system
- Selector switch for clean/city water
- Oversize junction box w/ master valves, filters & regulators
- Control mounted touch pad

Note: All Engle packages are standard with contour vinyl upholstery.